


ConcertoGEL Combs SELECTION GUIDE

ConcertoGEL Horizontal Electrophoresis Systems comprise a wide catalogue of casting combs for gel casting. Custom combs can be manufactured on request.



Combs for ConcertoGEL Series come in 4 different thicknesses and colour coded for ease of identification.

Black – 0.75 mm for tightly resolved bands White – 1 mm supplied as standard comb
 Red – 1.5 mm to maximise sample volume Blue – 2 mm to maximise sample volume

Black and white combs are recommended for high resolution gels and publication quality data. Red and blue are to scale-up nucleic acid volumes for preparatory techniques.


ConcertoGEL Midi




CGMD1-1, (360µl)




CGMD2-1, (158µl)




CGMD4-1, (77µl)




CGMD8-1, (41µl)




CGMD10MC-1, (27µl)




CGMD12-1, (23µl)



CGMD16-1, (16µl)



CGMD20-1, (14µl)



CGMD25-1, (10µl)

ConcertoGEL Classic



CGCL1-1, (495µl)



CGCL2-1, (225µl)




CGCL4-1, (122µl)




CGCL10-1, (45µl)




CGCL10MC-1, (29µl)



CGCL12-1, (41µl)



CGCL14MC-1, (29µl)



CGCL16MC-1, (27µl)



CGCL20MC-1, (10µl)




CGCL28MCA-1, (11µl)




CGCL35-1, (10µl)


ConcertoGEL Maxi




CGMX1-1, (675µl)




CGMX2-1, (315µl)




CGMX4-1, (153µl)




CGMX10-1, (72µl)




CGMX16-1, (41µl)




CGMX20MC-1, (27µl)




CGMX25-1, (21µl)




CGMX30-1, (17µl)



CGMX36-1, (14µl)



CGMX40MCA-1, (11µl)



CGMX50-1, (10µl)



MC = multichannel pipette compatible



MCA = fits multichannel loading every other row