

# Hercuvan

Lab Systems

## Disposable Gel Cutting Tips

### HIGHLIGHTS

- Specially designed to excise bands from agarose and protein gels
- Easy and convenient with one-hand operation
- Avoid scalpel damage to gel tray and transilluminators
- Push button gel and tip release mechanism
- Suitable for standard 1000µl pipettes
- Autoclavable, sterile versions also available
- Certified DNase-free, RNase-free and Pyrogen-free



## Convenient and safe for gel cutting

Hercuvan Disposable Gel Cutting Tips are specially designed to offer convenience and user-friendly operation to excise nucleic acid and protein bands from gels. It allows safe band removal and free from contamination. Furthermore, using sharp scalpels may damage the transilluminators and gel trays. The gel band is neatly extracted and contained in the tip, followed by dispensing into a tube for downstream application or analysis.

### ORDERING INFORMATION

Cat. No.	Disposable Gel Cutting Tips
LW-GCT4	Disposable Gel Cutting Tips, 4mmx 1mm, 250pcs/bag
LW-GCT4RACK	Disposable Gel Cutting Tips, 4mmx 1mm, 5 racks of 48pcs
LW-GCT6.5	Disposable Gel Cutting Tips, 6.5mmx 1mm, 250pcs/bag
LW-GCT6.5RACK	Disposable Gel Cutting Tips, 6.5mmx 1mm, 5 racks of 48pcs