

ElitePro FlexiStand

HIGHLIGHTS

- Fits wide range of pipettes
- “Windows” to accommodate multi-channel pipettes with ease
- Non-skid rubber bottom for stability
- Robust build with 5mm acrylic (40% thicker)

Versatile and trusty companion for your pipettes

ElitePro FlexiStand is designed to fit both single and multi-channel pipettes with reliable stability. The stand features “windows” to accommodate multi-channel pipettes head longways through the rack while simultaneously keeping the body of the pipette stable, which is achieved with two points; under the finger hook and near the bottom of the handle. As a result, users do not have to worry about pipettes falling or damaged. ElitePro FlexiStand is available as a 6-place rack that accepts four multi-channel pipettes. Each ElitePro FlexiStand is fitted with four non-skid rubber bottom for extra stability.



TT-EP-FS6

ElitePro FlexiStand features windows which can accommodate multi-channel pipettes

ORDERING INFORMATION

Cat. No.

ElitePro FlexiStand Rack

TT-EP-FS6

ElitePro FlexiStand for 6 pipettes, up to 4 multi-channels