

X-Blu UHC Western Blot Film

HIGHLIGHTS

- High performance, contrast and sensitivity
- Low background
- Light blue base
- Publication quality
- Suitable for ^{14}C , ^{32}P , ^{33}P , ^{125}I , ^{35}S

High Performance, Contrast and Sensitivity

X-Blu UHC Western Blot Film is an autoradiography film suitable for detection and documentation of ECL-based Western Blots, chemiluminescent nucleic acid or protein applications such as Northern and Southern blotting. The light blue base X-ray film provides optimized balance of ultra-high contrast and sensitivity with low background. It is also coated on both sides (double emulsion) which provides enhanced quality and performance. Clear and sharp images of publication quality can be obtained with X-Blu UHC Western Blot Film.



ORDERING INFORMATION

Cat. No.	X-Blu UHC Western Blot Film
LW-UHC5x7	X-Blu UHC Western Blot Film, 5 x 7 in. (127 x 178 mm), 100 sheets
LW-UHC8x10	X-Blu UHC Western Blot Film, 8 x 10 in. (203 x 254 mm), 100 sheets